

# Left-Hand Drive Syringe Pump System

## NE-M500 Complex Simplified™ by New Era Pump Systems Inc.

A simple, cost-effective - and versatile, pumping system to position a syringe to dispense towards the **left, rear, up, down - wherever**, while maintaining a **forward-facing keypad and display**.

For constrained work spaces that require positioning of a syringe pump while maintaining access to the keypad and display.

### The NE-M500 Base Package Includes:

- NE-500 OEM Single channel syringe pump
- Bottom cover assembly
- NE-1100X Programmable control box with upgraded 'X' version firmware
- CBL-DUAL-KIT containing all necessary cables to connect modules to each other and to a computer
- Power supply
- PWRSUP-EXT-2 Power supply 2 plug extension cable
- User manuals

# \$1150



The control box can be programmed and setup as any of our programmable syringe pumps, including the uploading of programs and remote control.

The base package can be **upgraded** to include **any OEM syringe pump model** and/or add additional pumps in a pump network configuration. For example, **the NE-500 one channel syringe pump** can be substituted with the NE-4500 two channel syringe pump. Or, both the NE-500 and NE-4500 can be networked together and operated as a 3-syringe system, operated by the control box.

## Place **Pump** in Any Direction!



**Smart • Simple • Cost-Effective**



**NEW ERA**  
Pump Systems Inc.  
SyringePump.com

138 Toledo Street  
Farmingdale,  
New York 11735



+1.631.249.1392



SyringePump.com



Info@SyringePump.com