


# Model: P-DKIT

Dual Pump Plumbing Kit

Functional Description

[www.SyringePump.com](http://www.SyringePump.com)



## P-DKIT Includes:

- (2) Dual Check valves (Part #Q-DVALVE)
- (2) 50 mL Air-Tite® syringes (Part #SYR-AT50LL)
- (2) "Y" Connectors (Part #Q-Y)
- (2) Female Luer to barb connectors (Part #Q-FL)
- (2) Male Luer to barb connectors (Part #Q-ML)
- Assembled with 0.125" ID PVC tubing (Part #Q-T125):  
6" connects the valves to the "Y" connectors, 5' to connect to your reservoir and out flow

- Withdraw pulls fluid from a reservoir
- Infusion out flows to a separate line
- Comes pre-assembled: Just attach the syringes
- Parts can be purchased individually

NE-1000 syringe pumps and synchronization cable CBL-DUAL-3 available separately.