



New Era
instruments
By New Era Pump Systems

Fast, Easy, Precise.

SyringeONE by New Era Instruments



+1.631.249.1392



info@SyringePump.com

The **SyringeONE*** is a single channel syringe pump capable of infusing or withdrawing fluids at extremely low or high rates. It contains features not available in the **InfusionONE**, such as stall detection, programmable dispense volumes, and more.



Applications



4-Step Setup

1. Set the syringe
2. Set the pumping
3. Select the direction
4. Insert the syringe and

Highlights

- **Simple interface** to learn and operate.
- Pumping Precision, within +/- 1 % accuracy.
- **Infusion and withdrawal** rates from **.73 µL/hr to 2120 mL/hr**.
- Infuses and withdraws fluids with **celerity and reliability**.
- Supports syringes from 1 mL to 140 mL.
- **Easy to program** simple dispenses, continuous flow, complex automation and more via the keypad or a PC.
- Available in a **High Pressure** version, which has double the force of the standard **SyringeONE**.
- Space saving footprint; only 8.75" x 5.75".
- **2 year warranty** and **lifetime technical support**.
- Can be **upgraded** with advanced firmware versions.

* This product is not an FDA cleared device.

Looking for more pumping features?

Visit www.SyringePump.com to view our wide selection of syringe pumps and the array of features and functions they offer.

MADE IN
USA